



SaniSure™ Accessory Kit-STR

Stirred Tank Reactor

Engineered Single Use Solutions...Reducing Operational Cost, Improving Production Efficiency



- Reduce Set-Up Time & Costs
- Consistency in Design Reduces Operator Error
- Easy Identification – All Lines Labeled
- Complete Assembled Kits or Individual Components Available
- Let SaniSure Customize To Your Needs

Pre-Fabricated Single-Use Bioreactor Assembly Kits

<u>Features</u>	<u>Benefits</u>
◇ Single-Use	<ul style="list-style-type: none"> ○ Minimize labor costs ○ Eliminate cleaning validation ○ Time installation savings - Reducing set up time and costs
◆ Pre-Fabricated Assemblies	<ul style="list-style-type: none"> ● Minimize injuries associated with building/cutting in-house, improving safety ● Improved consistency by eliminating operator dependency
◇ Offers Closed System	<ul style="list-style-type: none"> ○ Gives ability to assemble in open lab type environment ○ Provides aseptic solution to close off an unused line, making it usable again
◆ Improved Performance	<ul style="list-style-type: none"> ● Minimize contamination using advance components and latest connective technology
◇ Full Documentation Package	<ul style="list-style-type: none"> ○ Kits fabricated with full traceability and accompanying validation documentation ○ Meets UPS Class VI requirements
◆ Tubing Material Options	<ul style="list-style-type: none"> ● Silicone, C-Flex, PharMed, or PVC materials based on individual application needs ● Provides individual different options for welding, pumping, clarity, etc.




SaniSure™
Be sure.

www.SaniSure.com

SaniSure Inc.
P.O. Box 2360
Moorpark, CA 93021
Toll Free (877) 777-3005

SaniSure is
a Trademark
of SaniSure Inc.



SaniSure Part Number	Description
AK-STR-COMPLETE	Complete Kit; 1 of Each Assembly- 3 of the Tri-port Line Assembly
AK-STR-MEDIA	Media Inlet Line Assembly
AK-STR-SPARG	Sparger Line Assembly
AK-STR-COND	Condenser/Vent Line Assembly
AK-STR-SAMPLE	Sample Line Assembly
AK-STR-TRIPORT	Tri-port Assembly - 3 Assemblies Included
AK-STR-OVERLAY	Overlay Assembly, Gas Inlet to Head Space - 1 Assembly Included
AK-STR-TPSING	Single Tri-port Replacement Assembly - 1 Assembly Included

Research Applications

- Process optimization
- Screening studies
- Product optimization
- Scale up/scale down studies
- Educational programs



Microbial & Cell Culture Suites

- Yeast Fermentation
- Bacterial Fermentation
- Cell Suspensions
- Mammalian and Plant Cell Culture

Sterilization:

- Autoclave Max – 15 psi @ 134°C, 60 min
- Gamma Radiation – 25 kGy to 50 kGy



Represented by:

For More Information Visit: www.SaniSure.com

IMPORTANT: It is the user's responsibility to ensure the suitability and safety of SaniSure, Inc. products for all intended uses.

For a period of 90 days from the date of first sale, SaniSure, Inc. warrants this product to be free from defects in materials and workmanship. Our only obligation will be to replace any portion proving defective, or at our option, to refund the purchase price thereof. User assumes all other risk, if any, including the risk of injury, loss or damage, direct or consequential, arising out of the use, misuse, or inability to use, this product. THIS WARRANTY IS IN LIEU OF THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE, AND ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. No deviation is authorized.

SaniSure, Inc. assumes no obligations or liability for any advice furnished by it or for results obtained with respect to those products. All such advice is given and accepted at the buyer's risk. SaniSure™ is a trademark of SaniSure, Inc. ©2009 All rights reserved